OPAK

PORTUGUÊS

INTRODUÇÃO

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ENGLISH

INTRODUCTION

With the aesthetic needs improvement on Dentistry, opakers composites were developed to cover metal and dark surfaces.

DEFINITION

OPAK is a light-cured masking resin for the covering of dark dentin, metal and carbon fiber posts. And for the covering of metal surfaces in the dental laboratory.

INDICATIONS

- Dark and stained teeth;
- Metal surfaces during repair of metal-ceramic crowns and other restorations with a metal framework;
- Metal posts and cores in teeth to be restored with metal-free ceramics, indirect composite resin or ceromer;
- Carbon fiber posts;
- Metal surfaces in the dental laboratory.

CONTRAINDICATION

OPAK is contraindicated in patients with history of allergic reaction to any of its components.

COMPOSITION

Bis-GMA, urethane dimethacrylate, catalysts, stabilizers, pigments.

COLORS

B 0,5*	A 3*	Light pink	Dark pink
--------	------	------------	-----------

** B and A are basic Vita shades.Vita is a trademark of Vita Zahnfabrik H. Rauter GmbH & Co.*

PROPERTIES

- High masking capacity;
- Low viscosity: can be applied with brushes;
- Excellent adhesion;
- Compatible with all bonding agents and composite resins.

INSTRUCTIONS FOR USE IN THE DENTAL OFFICE

E. MASKING OF DARK DENTIN

- Isolate the operative field (pic. 1);
- Proceed with:
 - Partial removal of dark dentin;
 - Application of ANGELUS® 37% PHOSPHORIC ACID for 15 seconds; wash and gently dry;
 - Application of bonding agent on dentin and light-curing.
- Select an OPAK shade (B 0.5 or A3) and paint (with spatula or brush) a thin layer (0.3 mm) on dark dentin (pic. 2);
- Light-cure for 40 seconds (pic. 3);
- If necessary, apply a second layer of OPAK for a more intense masking effect or to change the initial masking shade (example: A3 layered on top of B 0.5). Repeat light-curing (pic. 4);
- Apply a composite resin and finish the restoration (pic. 5).

F. MASKING OF METAL POSTS AND CORES

- Isolate the operative field (pic. 1);
- Increase retentivity of metal (roughening with burs or sandblasting) or use a primer for metal (MZ PRIMER ANGELUS®) (pic. 2);
- Apply a bonding agent and light-cure (pic. 3);
- Select an OPAK shade (B 0.5 or A3) and paint (with spatula or brush) a thin layer (0.3 mm) on metal (pic. 4);
- Light-cure for 40 seconds;
- If necessary, apply a second layer of OPAK for a more intense masking effect or to change the initial masking shade (example: A3 layered on top of B 0.5). Repeat light-curing (pic. 5);
- Final result (pic. 6).

G. MASKING OF METAL IN CROWN REPAIRS

- Isolate the operative field; bevel the fractured ceramic or indirect resin surface and roughen the exposed metal with a diamond bur (pic. 1);
- Apply 10% ANGELUS® PORCELAIN ETCHANT gel of hydrofluoridric acid (HF) on the beveled surface and on the fracture site of ceramics or composite for 1 minute; wash abundantly with water and dry (pic. 2);

❶ IMPORTANT: Protect oral tissues against the corrosive and poisonous effects of HF.

- Apply SILANE ANGELUS® with brush on ceramics or composite; wait for 1 minute and gently dry with air (pic. 3a);
- Apply a bonding agent and light-cure (pic. 3b);
- Select an OPAK shade (B 0.5 or A3) and paint (with spatula or brush) a thin layer (0.3 mm) on metal only. Light-cure for 40 seconds (pic. 4);
- If necessary, apply a second layer of OPAK for a more intense masking effect or to change the initial masking shade (example: A3 layered on top of B 0.5). Repeat light-curing (pic. 5);
- Apply a composite resin and finish the repair (pic. 6).

H. MASKING OF CARBON FIBER POSTS

- After field isolation, dry the cemented post with air (pic. 1);
- Select an OPAK shade and paint the post with a thin layer (0.3 mm) (pic. 2);
- Light-cure all surfaces for 40 seconds each. If a more intense masking effect is needed, apply a second layer of OPAK. Repeat light-curing (pic. 3);
- Composite resin core build-up ready (pic. 4).

INSTRUCTIONS FOR USE IN THE DENTAL LAB

- Sandblast and apply a metal primer on the metal surface MZ PRIMER ANGELUS® (pics. 1 and 2);
- Apply a thin layer of the selected OPAK (light or dark) with a brush (pic. 3);
- Light-cure for 3 minutes;
- Apply two extra layers; light-cure each for 3 minutes (pic. 4);
- Light-cure the final layer for 6 minutes (pic. 5).

WARNINGS

- Replace syringe cap immediately after use. Prolonged exposure to any source of light may start polymerization;
- Contact with non-polymerized OPAK may irritate oral tissues, eyes and skin and cause sensitization in susceptible individuals;
- Keep OPAK in temperatures between 2°C and 28°C.

PRECAUTIONS

- Do not apply thick layers of OPAK;
- Isolate area ideally with rubber dam;
- Avoid contact of OPAK with eugenol-based substances which inhibit its polymerization;
- Prepare the metal surface before the application of OPAK (sandblasting and metal primer - MZ PRIMER ANGELUS®);
- Keep OPAK in refrigerated environment, with temperature lower than 28°C;
- Replace syringe cap immediately after use. Prolonged exposure to any source of light may start polymerization;
- Contact with non-polymerized OPAK may irritate oral tissues, eyes and skin and cause sensitization in susceptible individuals.

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